


# New!!! NEAR ZERO CARRY OVER

Meeting today's Mass Spectrometers sensitivity requirements

- Fast** *Typical Clean Cycle time less than 1 minute*
- Near Zero Carryover** *Typically less than 0.003% (30ppm)*
- Minimum Sample Contact** *Wetted parts are SS & PEEK*
- Reproducible** *Better than 1% RSD with partial loop fill*
- Wash** *the retention tubing, needle and valve all at once with up to 2 solvents*
- Upgradable** *DLW Upgrade kit available for existing PALs*

## Carryover Results UV

- 2ul Injection Chlorhexidine** “notoriously sticky substance”
- 1200ng Chlorhexidine on column, followed by 2 blank Injections**
- Detection UV 257nm**



Conditions for UV carryover test

Column:	Halo C18 2.1x50mm, 2.8 $\mu$ m
Flow:	0.5ml/min
Injection:	Full Loop 2 $\mu$ l (PEEK)
Eluent:	A: H <sub>2</sub> O + 0.1% TFA / B: Acetonitrile + 0.1% TFA
Gradient:	10%B (0.5 min) – 10% to 90% B (7.5 min) – 90% B (1 min)
Wash1:	H <sub>2</sub> O + 0.1% TFA
Wash2:	Acetonitrile + 0.1% TFA
Detection:	UV, 257nm

