



## ITEX-2

### Application Type

Standard Product

### Application ID

Extend the headspace mode

### Description

Accessory Showcase

#### Contents [\[hide\]](#)

- [1 ITEX-2 Option](#)
- [2 Significant Markets](#)
- [3 Steps of CombiPAL ITEX Cycle](#)
- [4 Photos](#)
- [5 Downloads](#)
- [6 PRINT Function for Page](#)
- [7 Accessories for the PAL](#)
- [8 Contact LEAP](#)



## ITEX-2 Option [\[edit\]](#)

The ITEX-2 Option consists of an add-on module which can be used with any existing or new CombiPAL robot. It performs enrichment of volatile or semi-volatile compounds during headspace analysis. A microtrap filled with adsorbent material, such as Tenax or activated charcoal is placed between the heated CombiPAL Headspace syringe and the syringe needle. No cryofocussing is needed to obtain sharp peaks. The ITEX-2 Option offers faster desorption, a complete deactivated sample path and a shortened syringe. Adding ITEX to your existing COMBI PAL autosampler is now as easy as switching from Liquid to Headspace Sampling.

- Rapid and efficient sample enrichment of volatile and semi-volatile compounds in solid, liquid and gaseous samples
- In-tube extraction and direct thermal desorption using proven industry standard adsorbents (Tenax TA, Carbotrap, Carbosieve etc.)
- Syringe only concept for transparent sample handling, no sample loops, transfer lines, or switching valves
- Top mounted on GC's, saves valuable bench space
- No GC injector modifications, no cryo-focussing necessary
- Interfaces with any CombiPAL System, controlled by all major GC/GC-MS

## Significant Markets


[edit]

- Environmental
- Drinking Water
- Foods / Flavors
- Consumer Products
- Forensics / Toxicology
- Petrochemicals
- Polymers
- Pharmaceuticals
- Residual Solvents

## Steps of CombiPAL ITEX Cycle

[edit]

- 1 Sample is heated and / or agitated in a sealed sample vial until equilibrium is achieved
- 2 The Trap needle pierces the vial septa and the syringe pumps the headspace through the Trap
- 3 The loaded ITEX Trap is flash heated up to 350°C and desorbed into the hot GC Injector
- 4 After thermal desorption the hot ITEX Trap is cleaned with inert flush gas



**Extraction Trap:** 44 mg adsorbent Tenax TA 80/100 mesh

**Pumping syringe:** 1.3ml HD syringe with removable trap

**Extraction speed:** 10µl/sec. - 1000µl/sec

**Extraction strokes:** 1 -999

**Extraction volume:** 130µl - 1300 µl per stroke

**Desorption temp:** +5°C above ambient – 350°C selectable in 1°C steps

**Desorption time:** 0 - 300 sec.

**Desorption speed:** 1µl/s - 500µl/s

**Heating-up rate:** up to 12°C/s

**Trap cleaning:** Inert gas purging, 30 - 600 sec.

**Sample capacity:** up to 297 2ml standard vials up to 96 10ml or 20ml vials optional Cooled Trayholders available

**Incubator oven:** 6 heated positions for 2ml/10ml/20ml vials


**Incubation temp:** 30°C – 200°C selectable in 1°C steps

**Incubation time:** up to 999 min.

**Agitation:** Interval shaking 250rpm – 750 rpm prior sample extraction

## Photos

[\[edit\]](#)



## Downloads

[\[edit\]](#)



[ITEX](#)



[ITEX Application Binder](#)

## PRINT Function for Page

[\[edit\]](#)



[PDF of this page formatted for printing](#)

## Accessories for the PAL



[Other Accessories for PAL Robots](#)

**For the latest up to date information and more photos and videos please visit:**

<http://www.leapwiki.com>